



**USAID**  
FROM THE AMERICAN PEOPLE

# Chest X-ray Taking Procedures Training for X-ray technicians/ Radiographer

## “X-ray Production”

Ms. KHIN YADANAR MOE

Consultant (TB CXR Training), IDDS Project/ Myanmar

# Content



Introduction



Equipment Used in Production of X-rays



Production of X-rays



X-ray Beam Characteristics



Quality Control of X-ray Producing Equipment



X-ray Tube Faults

# Introduction

“X-rays are used in medical imaging to produce the images by penetrating the **high energy** x-rays photons into the internal structures of the body and captured on the image receptor (Analog and Digital acquisition).”

Reference: X-ray technician/radiographer TB Chest X-ray Training Curriculum by Prof U Khin Hla

# Introduction – Cont.

After passing the X radiation through the body:


Bone (dense) - white  
Fat - dark  
Soft tissue - grey  
Air - black

Soft tissue

Bone

Fat

## CXR PA View



Reference: X-ray technician/radiographer TB Chest X-ray Training Curriculum by Prof U Khin Hla

# Equipment Used in the Production of X-rays

1. Autotransformer (A variable transformer)
2. Voltmeter
3. Timer
4. X-ray tube
5. Focal spot
6. Exposure switch
7. Glass envelope
8. Tube housing

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.

# Equipment Used in the Production of X-rays – Cont.

**Autotransformer  
(A variable  
transformer)**



**Voltmeter**



**Timer**



Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.

# Equipment Used in the Production of X-rays – Cont.

## Electronic timer

- Contained in most radiographic equipment
- Allows exposure times of 1 ms (0.001 second)

## Automatic exposure control (AEC)

- To terminate the exposure when a predetermined amount of exposure has been reached
- To provide consistent exposure (film, digital detector)
- Shortest time with an AEC is 1 ms (0.001 second)
- AEC sensor is placed between the patient and the image receptor

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong ScD FACR FACMP, Stewart.**

# Equipment Used in the Production of X-rays – Cont.

## X-ray tube



## Focal spot



Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.



# Equipment Used in the Production of X-rays – Cont.

## Exposure switch



## Glass envelope



## Tube housing



Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.

# Production of X-rays

Three types of conditions:

- i. Source of electrons
- ii. Acceleration of electrons
- iii. Sudden stoppage of electrons against target material

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.

# Production of X-rays – Cont.

*Note: If the machine has been **off overnight**; **warm-up exposures are needed** to warm the anode throughout (anode cracking can occur when surface heat is applied to a cold anode)*



Activating the filament circuit and heating the x-ray tube filament



To warm the filament and ready it for much higher current



Use technique charts for different examination

*Note: Machines having a line-voltage compensator on the control panel should be adjusted to **compensate for** any incoming voltage **fluctuation***

# Production of X-rays – Cont.

**two-stage exposure button**

It is the induction motor bringing anode rotation up to speed

Then the filament is heated to maximum (thermionic emission) and produces an electron cloud

The **first click** heard after partial depression

It leads to the exposure being made


The **second click** - after complete depression of the rotor/exposure switch

- Voltage selected by the autotransformer is sent to the step-up transformer
- Converted to the high voltage (kV)
- Low amperage (mA) required

- High-voltage current then passes through the rectification system that changes AC to pulsating DC

- Applied high voltage (potential difference) propels the electron cloud to the anode

# Production of X-rays – Cont.




# Types of Technique Charts

(a) Fixed kVp–variable mAs	(b) Variable kVp–fixed mAs	(c) Variable technique (vary both mAs and kVp)
Assumes optimal kVp for the part being radiographed	kVp is varied according to part thickness	Provides for alteration of routine techniques because of pathologic conditions, patient age, body mass index, contrast media
mAs is varied according to the part thickness		

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.

# Control Console



# Production of X-rays – Cont.

## Brems


- **X-rays** produced by slowing of incoming electrons by the target atoms; slowing releases energy in the form of x-rays

## Characteristic

- **X-rays** produced when incoming electrons at the anode dislodge orbital electrons from the target material, and outer shell electrons fall in to fill the hole created; this movement releases energy in the form of x-rays



# Production of X-rays – Cont.



## Production of X-ray beam

Reference: X-ray technician/radiographer TB Chest X-ray Training Curriculum by Prof U Khin Hla

# X-ray Beam Characteristics

- The resultant x-ray beam contains many different energies and is heterogeneous
- The maximum energy an x-ray photon can have corresponds to the kVp that was used

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.

# X-ray Beam Characteristics – Cont.

- Beam characteristics may be **altered** by using **filtration**:

## Inherent filtration

- the oil and glass envelope of the x-ray tube

## Total beam filtration

- equals **inherent filtration plus added filtration** (at least 2.5-mm aluminium equivalent)

Advantages:  
Filtration removes the low-energy (long-wavelength) rays from the beam and **reduces the patient dose**

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong ScD FACR FACMP**, Stewart.

# X-ray Beam Characteristics – Cont.

## Aluminum filter

A filter is usually a **sheet of aluminium** placed in the primary beam just as it exits the x-ray tube and before it reaches the collimator

## Other types of filters

Compensating filters (e.g., wedge, boomerang)

Half-value layer:


Amount of filtration that reduces the beam intensity by half

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong ScD FACR FACMP**, Stewart.

Low energy x-rays

Filter

High energy x-ray



## Filtration

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong ScD FACR FACMP**, Stewart.

# Quality Control of X-Ray–Producing Equipment

To provide safe and reliable operation of equipment:

- *Filtration-beam quality:*

Tested using a digital dosimeter

Half-value layer measurement is required

- *Collimator/light field to radiation field alignment:*

Must be accurate within 2% of SID (0.8” at 40” SID)

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong ScD FACR FACMP, Stewart.**

# Quality Control of X-Ray–Producing Equipment – Cont.

- *Effective focal-spot size:*
  - Should be within 50% of size stated in equipment specifications
- *kVp:*
  - Must be accurate to within 10% of that chosen
- *Timer:*
  - Should be within 5% of the time chosen for exposures over 10 milliseconds

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong** ScD FACR FACMP, Stewart.

# X-ray Tube Faults

No	Faults	Symptoms	Possible Causes
1	Reduced insulation oil	- An arcing sound heard during an exposure	- Repeated over heating or leakage of oil (forming air bubbles)
2	Deposition of vaporized tungsten on the glass wall	- Glass wall becoming colored based on age of use	- Heavy exposure increases vaporization of tungsten

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong ScD FACR FACMP**, Stewart.



# X-ray Tube Faults – Cont.

<b>No</b>	<b>Faults</b>	<b>Symptoms</b>	<b>Possible Causes</b>
3	Glass puncture	- Oil entering the tube	- Electric discharge through oil at high kv when oil insulation is inadequate
4	Anode wobbling	- Apparent movement of light patch on fluoroscopic screen	- Anode stem being bent due to excessive heat conduction during prolonged fluoroscopy

Reference: 10e **RADIOLOGIC SCIENCE FOR TECHNOLOGISTS: PHYS, BIOL PROTECTION** By **Bushong ScD FACR FACMP**, Stewart.

**THANK YOU !**