


# Mycobacteria Product Series: GenoType Test Systems

## (A) DNA Extraction Area, BSL 2/BSL 3 Laboratory

- Put on new gloves and lab coat.
- Decontaminate work area with freshly diluted 1.5% sodium hypochlorite solution.
- If Internal Control DNA (IC) is provided with the **GenoType** kit (**GenoType Mycobacterium CM VER 2.0/AS**, **GenoType CMdirect**), then prepare A-LYS/IC mix (A-LYS provided with **GenoLyse**® kit) for all samples and negative control first. Use 100 µl Lysis Buffer (A-LYS) and 2 µl IC per sample.
- If a negative control shall be included, use 100 µl A-LYS or A-LYS/IC mix and continue with **GenoLyse**® procedure, step 3.

### Starting material

#### GenoType MTBDRplus/sl VER 2.0 GenoType CMdirect



Patient specimen:  
Use 500 µl  
decontaminated  
patient specimen.



#### GenoType MTBDRplus/sl VER 2.0, GenoType Mycobacterium CM VER 2.0/AS, GenoType MTBC, GenoType NTM-DR




Bacteria grown in  
liquid medium:  
Use 1 ml  
of culture.



Bacteria grown on solid  
medium: Collect with  
inoculation loop, suspend in  
100 µl A-LYS or A-LYS/IC mix,  
vortex and continue with  
**GenoLyse**® procedure, step 3.



### DNA extraction with GenoLyse®



1. Centrifuge for 15 min at 10,000 x g in a centrifuge with aerosol-tight rotor, carefully remove supernatant.
2. Add 100 µl A-LYS or A-LYS/IC mix to pellet and vortex.
3. Incubate for 5 min at 95°C in a water bath, briefly spin down.
4. Add 100 µl Neutralization Buffer (A-NB) and vortex for 5 sec.
5. Centrifuge for 5 min at full speed.
6. = DNA solution  
For long-term storage, transfer supernatant to a new tube and store at -20°C.

- Decontaminate pipettes, tip boxes, racks and work area with freshly diluted 1.5% sodium hypochlorite solution.


Please make sure that the procedure described in this flow chart corresponds to the instruction for use (IFU) of the kit at hand. To do so, go to [www.hain-lifescience.de/ifu.html](http://www.hain-lifescience.de/ifu.html) and enter the IFU number of the test you are using.

# Mycobacteria Product Series: GenoType Test Systems

## (B) Pre-PCR Area, BSL 1 Laboratory


- Put on new gloves and lab coat.
- Decontaminate work area with freshly diluted 1.5% sodium hypochlorite solution.
- Avoid handling and pipetting above open vials, close lids after each pipetting step.
- Whenever opening new AM-A or AM-B, aliquot and store at -20°C.

### Prepare master mix containing AM-A and AM-B (per sample)



- Master mix needs to be prepared freshly each time.
- Do not vortex master mix.
- Determine the number of samples to be analysed plus controls and prepare master mix accordingly.

### Aliquot master mix into PCR tubes



- Decontaminate pipettes, tip boxes, racks and work area with freshly diluted 1.5% sodium hypochlorite solution.


Please make sure that the procedure described in this flow chart corresponds to the instruction for use (IFU) of the kit at hand. To do so, go to [www.hain-lifescience.de/ifu.html](http://www.hain-lifescience.de/ifu.html) and enter the IFU number of the test you are using.

# Mycobacteria Product Series: GenoType Test Systems

## (C) PCR Area, BSL 1 Laboratory


- Put on new gloves and lab coat.
- Decontaminate work area with freshly diluted 1.5% sodium hypochlorite solution.

### Add DNA solution to each master mix aliquot



- Volumes described are applicable for DNA extracted with **GenoLyse®**.
- For negative control use 5 µl water or 5 µl negative control isolate instead of DNA solution.
- If a positive control shall be included, pipette 5 µl C+ (**GenoType Mycobacterium CM VER 2.0/AS**, **GenoType CMdirect**) to one aliquot. Open and pipette C+ last and close vial immediately afterwards.

### Place PCR tubes into thermal cycler



- When using a thermal cycler with the respective pre-installation, select specific protocol for patient specimens or cultivated samples.

- Decontaminate pipettes, tip boxes, racks and work area with freshly diluted 1.5% sodium hypochlorite solution.

Please make sure that the procedure described in this flow chart corresponds to the instruction for use (IFU) of the kit at hand. To do so, go to [www.hain-lifescience.de/ifu.html](http://www.hain-lifescience.de/ifu.html) and enter the IFU number of the test you are using.


# Mycobacteria Product Series: GenoType Test Systems

## (D-E) Hybridization and Detection Area, BSL 1 Laboratory

- Put on new gloves and lab coat.
- Decontaminate work area with freshly diluted 1.5% sodium hypochlorite solution.
- Pre-warm RIN and distilled water to room temperature (RT) and HYB and STR to 45°C.
- Freshly dilute CON-C 1:100 with CON-D and SUB-C 1:100 with SUB-D in the amounts needed.



Use **TwinCubator** with program for the hybridization and detection procedure



- Decontaminate pipettes, tip boxes, racks and work area with freshly diluted 1.5% sodium hypochlorite solution.
- Please make sure that the procedure described in this flow chart corresponds to the instruction for use (IFU) of the kit at hand.  
To do so, go to [www.hain-lifescience.de/ifu.html](http://www.hain-lifescience.de/ifu.html) and enter the IFU number of the test you are using.